

Technical Data Sheet

TAE BUFFER 50X

Product Name: **TAE Buffer 50X Sterile Stock Solution**
Product Code: **TAExxxx**
Storage Conditions: **Room Temperature, Store Away From Light.**
Shelf Life: **2 years from date of manufacture.**

Content

Tris (2M)	242 g/L
Glacial Acetic Acid (1M)	60 g/L
Na ₂ EDTA (0.05M)	30.3 g/L

*Dilute 20ml TAE Buffer 50X with 980ml deionised water to make 1 litre of TAE Buffer.
Final concentration is 40mM Tris Base, 20mM Glacial Acetic Acid and 1mM Na₂EDTA.*

For further enquiries contact Formedium Technical Services.

Formedium Ltd products are supplied for Research Purposes only.