

## Technical Data Sheet

---

### TBE BUFFER 10X

Product Name: **TBE Buffer 10X Sterile Stock Solution**  
Product Code: **TAExxxx**  
Storage Conditions: **Room Temperature, Store Away From Light**  
Shelf Life: **15 months from date of manufacture.**

#### Content

Tris	(0.9M)	108 g/L
Borate	(0.9M)	55 g/L
Na <sub>2</sub> EDTA	(0.02M)	43.5 g/L

*Dilute 100ml TBE Buffer 10X with 900ml deionised water to make 1 litre of TBE Buffer. Final concentration is 89mM Tris Base, 89mM Borate and 2mM Na<sub>2</sub>EDTA.*

*For further enquiries contact Formedium Technical Services.*

*Formedium Ltd products are supplied for Research Purposes only.*